

NEST Compatibilty Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μ L)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016		
Axygen	Axyjet	Single	0.1-2	•			•			•				
			0.5-10	•			•			•				
			2-20		•			•				•		
			5-50		•			•						
			10-100		•			•						
			20-200		•			•						
			100-1000				•			•			•	
		Multiple	1-10	•				•			•			
			5-50		•			•						
			20-200		•			•						
			50-300		•			•						
		Biohit	Proline	Single	0.1-2.5	•			•			•		
					0.2-10	•			•			•		
					0.5-10	•			•			•		
5-50					•			•						
5-100					•			•						
10-250					•			•						
10-500							•			•				
20-200					•			•						
50-1000							•			•				
100-1000							•			•			•	
Multiple	0.2-10			•				•			•			
	5-100				•			•		•				

NEST Compatibilty Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μ L)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016		
			25-250		•			•						
			50-250		•			•						
			5-100		•				•					
	e-Line	Single	0.1-10	•				•			•			
			0.2-10	•				•			•			
			5-120		•				•					
			10-300		•				•					
			50-1000			•				•			•	
		Multiple	0.2-10	•					•			•		
			5-120		•				•					
			10-300		•				•					
	m-Line	Single	0.1-3	•				•			•			
			0.5-10	•				•			•			
			10-100		•				•					
			20-100		•				•					
			100-1000			•				•			•	
		Multiple	0.5-10	•					•			•		
			5-100		•				•					
			30-300		•				•					
	BrandTech	Transferpette	Single	0.1-1	•			•			•			
0.5-10				•			•			•				
2-20					•				•			•		
5-50					•				•					

NEST Compatibilty Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μL)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016		
		Single	10-100		•			•						
			20-200		•			•						
			25-250		•			•						
			100-1000			•				•			•	
		Multiple	0.5-10	•			•					•		
			2-20			•				•			•	
			2.5-25			•				•				
			5-50			•				•				
			10-100			•				•				
			20-200			•				•				
			30-300			•				•				
		Capp Aero		Single/Multiple	0.2-2	•			•				•	
					0.5-10	•			•				•	
5-50						•			•					
10-100						•			•					
25-200						•			•					
50-300						•				•				
100-1000								•			•			•
Labnet	Labpette Discovery or Biopette™ Plus	Single	0.1-2	•			•				•			
			0.5-10	•			•				•			
			2-20			•			•				•	
			10-100			•			•					
			20-200			•				•				

NEST Compatibilty Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μ L)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016		
			100-1000			•			•			•		
		Multiple	1-10	•			•			•				
			5-50		•			•						
			20-200		•			•						
			50-300		•			•						
	BioPette E	Single	0.5-10	•			•				•			
			2-20	•			•				•	•		
			10-200		•			•						
			100-1000				•			•			•	
		Multiple	0.5-10	•				•				•		
			2-20	•				•				•	•	
			10-200		•			•						
			100-1200				•				•			•
	Matrix	Impact/ Impact ²	Single/Mul tiple	2-125		•			•					
Nichryo/ Oxford	Benchmate	Single	0.5-10	•			•				•			
			2-20		•			•			•			
			10-50		•			•						
Rainin	EDP3	Single	E-MIC-10,0.5- 10	•			•				•			
			E2-25,2.5-25		•			•						
			E2-100,10-100		•			•						
			E2-250,25-250		•			•						

NEST Compatibilty Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μ L)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016		
Socorex			E2-1000,100-1000			•			•			•		
		Multiple	EP-M8-10,0.5-10	•			•			•				
			EP-M8-250,10-250		•			•						
	Latch Mode	Single	R2,0.1-2	•			•				•			
			R10,0.2-10	•			•				•			
			R20,2-20		•			•				•		
			R100,10-100		•			•						
			R200,20-200		•			•						
			R1000,100-1000				•			•				•
	Socorex	Callbra/ Accura	Single	0.1-2	•			•			•			
				0.2-2	•			•			•			
0.5-10				•			•			•				
1-10				•			•			•				
2-20					•			•				•		
10-100					•			•						
20-200					•			•						
100-1000							•			•				•
Multiple			1-10	•				•				•		
			5-50		•			•						
			10-100		•			•						

NEST Compatibilty Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μL)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016		
			20-200		•			•						
			50-200		•			•						
VWR	VWR and VWR Signature	Single	0.1-2	•			•			•				
			0.5-10	•			•			•				
			2-20		•			•				•		
			5-50		•			•						
			10-100		•			•						
			20-200		•			•						
			50-250		•			•						
		100-1000				•			•				•	
		Multiple	1-10	•				•				•		
			5-50		•			•						
			20-200		•			•						
			50-300		•			•						
Eppendorf	Eletronic	Single/Mul tiple	0.5-10	•			•			•				
			5-100		•			•						
			2-200		•			•						
			50-1000				•			•			•	
	Reference	Single/Mul tiple	0.1-2.5	•			•				•			
			0.1-10	•			•				•			
			0.5-10	•			•				•			
			2-20		•			•				•		
			10-100		•			•						

NEST Compatilby Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μ L)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016		
Thermo			20-200		•			•						
			50-200		•			•						
			100-1000			•				•			•	
	Research	Single	0.1-2.5	•				•			•			
			0.5-10	•				•			•			
			2-20		•				•				•	
			10-100		•				•					
			20-200		•				•					
			100-1000			•				•				•
		Multiple	0.5-10	•				•				•		
			5-100		•				•					
			10-100		•				•					
			50-100		•				•					
			20-300		•				•					
			30-300		•				•					
			50-1200			•					•			
	Research Pro	Single	20-300		•				•					
			50-1000			•				•				
		Multiple	0.5-10	•				•			•			
			20-300	•				•			•			
Research Plus	Single/Multiple	100-1000			•			•			•			
Thermo	Finnpipette	Single	0.2-2	•			•			•				

NEST Compatibilty Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μL)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016	
Scientific			0.3-3	•			•			•			
			0.2-10	•			•			•			
			0.5-5	•			•			•			
			0.5-10	•			•			•			
			1-10	•			•			•			
			5-40		•			•					
			5-50		•			•					
			10-100		•			•					
			40-200		•			•					
			30-300		•			•					
			100-1000			•				•			•
			200-1000			•				•			•
	Multiple		1-10	•				•			•		
			5-50		•			•					
			50-300		•			•					
	Finnpipette Digital	Single	40-200			•							
Finnpipette electronic	Multiple	1-10	•				•			•			
		30-300		•			•						
		100-1000			•				•			•	
Gilson Pipetman	Pipetman Classic	Single	P2,0.1-2	•			•			•			
			P2,0.2-2	•			•			•			
			P10,0.2-10	•			•			•			

NEST Compatibilty Chart of Filter tips

Brand	Pipette	Channel	Channel Volume (μL)	311001	312001/ 316001	313001	311012	316012/ 312012	313012	314016	310012	313006/ 313016	
			P20,2-20		•			•			•		
			P100,10-100		•			•					
			P200,20-200		•			•					
			P1000,100-1000			•			•			•	
	Pipetman Neo	Single	P2N,0.2-2	•			•				•		
			P10N,1-10	•			•				•		
			P20N,2-20		•			•				•	
			P100N,10-100		•			•					
			P100N,20-100		•			•					
			P200N,20-200		•			•					
		Multiple	P8×20N,2-20		•			•					
	Pipetman Ultra	Single	U200,20-200		•				•				
		Multiple	P200-M8,20-200		•				•				

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

