

 ALLIANCE  
BIO EXPERTISE

 **RIEGER**  
Industrievertretungen GmbH  
[www.rieger-iv.at](http://www.rieger-iv.at)

Wir sind für Sie da:

E-Mail: [office@rieger-iv.at](mailto:office@rieger-iv.at)


Telefon: +43(1) 728 00 52

Fax: +43(1) 728 69 16

Rustenschacher Allee 10

1020 Wien, Austria

**AIRWEL** ®  
Bio-collector



**ДИА•М**  
современная лаборатория

[www.dia-m.ru](http://www.dia-m.ru)  
заказ on-line



**AIRWEL®**

Increase your productivity  
quality  
traceability

**ADAPTABLE TO YOUR NEEDS**

PETRI DISH:  
90 mm  
55 mm


AIRFLOW:  
100L / min  
200L / min

**SPECIFIC OPTION**

COMPRESSED GAS SYSTEM  
for Controlled Environment

CONNECTED  
MICROBIOLOGY

**FULL TRACEABILITY**  
GMP / GLP – 21 CFR PART 11



**IMPORT DATA**  
Bluetooth bar code scan

**EXPORT DATA**

**EXCLUSIVE OPTION**

**IN SITU AIRFLOW CHECK**

Check your airflow whenever you like!

Don't wait the 1 year calibration services to check

Avoid wrong collection due to wrong airflow!!

**SAFETY MANAGEMENT**

Cascade password  
3 level:  
Operator – Manager – Service



**FULLY HERMETIC**  
ip65 dust and water protection

**WIRELESS BATTERY CHARGER**

**ANTIBACTERIAL PAINT**

**SERVICES**  
MANUFACTURER or DEALER **CALIBRATION**  
QUALIFICATION: **IQ.OQ** protocol



**SPECIFICATIONS**

**AIRWEL®**



COMPLIANCE	AIRWEL	AIRWEL+
IP65		✓
ISO14698 / ISO17025	✓	✓
USP797 / USP116	✓	✓
EU GMP / GLP	✓	✓
21 CFR PART 11	✓	✓
TRACEABILITY		
Bluetooth module	✓	✓
PC software for data management	✓	✓
Android App (PC / Smartphone / tablet)	✓	✓
CALIBRATION		
IN SITU AIRFLOW CHECK (Option)		✓
A.B.E or dealer calibration	✓	✓
BENEFITS		
Cascade password		✓
Wireless Induction charger		✓
Battery autonomy	More than 10 hours	
Type of dishes	90 mm or 55 mm	
Type of heads	Stainless steel or Disposable polyester	
Airflow	100L or 200L / min	
Transport case	Light	Rigid
Multi-mode aspirating cycles	Uninterrupted	Fragmented
Antigerm paint	✓	✓
DIMENSIONS		
W x H x D	130 x 330 x 150 mm	
Weight	1.6 kg	
Power supply / Electric power	120-230 V – 50-60HZ / 120W	
Made in EU	✓	✓




# ALLIANCE BIO EXPERTISE

Empowered from our people's skills and experiences, we, ALLIANCE BIO EXPERTISE, together with our customers and partners decided to manufacture Microbiology automated equipment around essential values:  
Customer Focus, Service, Loyalty, Innovations

**Ask your local partner  
for a personal demonstration!**

**000 «Диаэм»**

**Москва**

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ [sales@dia-m.ru](mailto:sales@dia-m.ru)

[www.dia-m.ru](http://www.dia-m.ru)

**С.-Петербург**

+7 (812) 372-6040  
[spb@dia-m.ru](mailto:spb@dia-m.ru)

**Новосибирск**

+7 (383) 328-0048  
[nsk@dia-m.ru](mailto:nsk@dia-m.ru)

**Воронеж**

+7 (473) 232-4412  
[vrn@dia-m.ru](mailto:vrn@dia-m.ru)

**Йошкар-Ола**

+7 (927) 880-3676  
[nba@dia-m.ru](mailto:nba@dia-m.ru)

**Красноярск**

+7 (923) 303-0152  
[krsk@dia-m.ru](mailto:krsk@dia-m.ru)

**Казань**

+7 (843) 210-2080  
[kazan@dia-m.ru](mailto:kazan@dia-m.ru)

**Ростов-на-Дону**

+7 (863) 303-5500  
[rnd@dia-m.ru](mailto:rnd@dia-m.ru)

**Екатеринбург**

+7 (912) 658-7606  
[ekb@dia-m.ru](mailto:ekb@dia-m.ru)

**Кемерово**

+7 (923) 158-6753  
[kemerovo@dia-m.ru](mailto:kemerovo@dia-m.ru)

**Армения**

+7 (094) 01-0173  
[armenia@dia-m.ru](mailto:armenia@dia-m.ru)

