

# Pipet Tip Compatibility

Pipettor Brand/Make	Volume Range	Volume/Biologix Tip	20-0010	21-0010	20-0200	21-0200	100-1000µl	20-1000	100-1000µl	21-1000	20-1250	100-1000µl	20-1251	100-1000µl	21-1250	22-0010	0.5-10µl	22-0011	0.5-10µl	23-0010	1-200µl	23-0011	1-200µl	22-0200	1-200µl	22-0201	1-200µl	23-0200	1-200µl	23-0201	100-1000µl	22-1000	100-1000µl	22-1001	100-1000µl	23-1000	100-1000µl	23-1001			
Axygen Axytip (single channel)	0.1-2µl	✓	✓										✓	✓	✓	✓																									
	0.5-10µl	✓	✓										✓	✓	✓	✓																									
	2-20µl			✓	✓																																				
	5-50µl			✓	✓																																				
	10-100µl			✓	✓																																				
	20-200µl			✓	✓																✓	✓	✓	✓																	
	50-250µl																																								
100-1000µl																																									
Axygen Axytip (multichannel)	1-10µl	✓	✓										✓	✓	✓	✓																									
	5-50µl			✓	✓																																				
	20-200µl			✓	✓																✓	✓	✓	✓																	
	50-300µl			✓	✓																																				
Axygen Axytip-16	0.2-2µl																																								
	0.5-10µl																																								
	5-50µl																																								
Biohit Proline (single channel)	0.2-10µl	✓	✓										✓	✓	✓	✓																									
	5-100µl			✓	✓																																				
	10-250µl			✓	✓																✓	✓	✓	✓																	
	10-500µl					✓	✓	✓	✓	✓																															
	50-1000µl					✓	✓	✓	✓	✓																															
Biohit Proline (multichannel)	0.2-10µl	✓	✓										✓	✓	✓	✓																									
	5-100µl			✓	✓																✓	✓	✓	✓																	
	25-250µl			✓	✓																✓	✓	✓	✓																	
Biohit e-Line (single channel)	1-10µl	✓	✓										✓	✓	✓	✓																									
	5-120µl			✓	✓																																				
	10-300µl			✓	✓																																				
	50-1000µl					✓	✓	✓	✓	✓																															
Biohit e-Line (multichannel)	0.2-10µl	✓	✓										✓	✓	✓	✓																									
	5-120µl			✓	✓																																				
	10-300µl			✓	✓																																				
Biohit m-Line (single channel)	0.5-10µl																																								
	2-20µl			✓	✓																																				
	10-100µl			✓	✓																																				

# Pipet Tip Compatibility

	0.5-10µl	0.5-10µl	1-200µl	1-200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	0.5-10µl	0.5-10µl	0.5-10µl	0.5-10µl	1-200µl	1-200µl	1-200µl	1-200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	Volume/Biohit Tip	Volume Range	Pipettor Brand/Make
	20-0010	21-0010	20-0200	21-0200	20-1000	21-1000	20-1250	20-1251	21-1250	22-0010	22-0011	23-0010	23-0011	22-0200	22-0201	23-0200	23-0201	22-1000	22-1001	23-1000	23-1001			
			✓	✓										✓	✓	✓	✓						20-200µl	
					✓	✓	✓	✓	✓										✓	✓	✓		100-1000µl	
✓	✓									✓	✓	✓	✓										0.5-10µl	Biohit m-Line (multichannel)
		✓	✓											✓	✓	✓	✓						5-100µl	
		✓	✓											✓	✓	✓	✓						30-300µl	
✓	✓									✓	✓	✓	✓										0.2-2µl	
✓	✓									✓	✓	✓	✓										0.5-10µl	
		✓	✓																				5-50µl	
		✓	✓																				10-100µl	Capp Aero (multichannel)
		✓	✓											✓	✓	✓	✓						25-200µl	
		✓	✓											✓	✓	✓	✓						50-300µl	
				✓	✓	✓	✓	✓											✓	✓	✓		100-1000µl	
																							0.2-2µl	
																							0.5-10µl	384 Format (multichannel)
																							5-50µl	
✓	✓									✓	✓	✓	✓										0.1-1µl	
✓	✓									✓	✓	✓	✓										0.5-10µl	
																							2-20µl	
		✓	✓																				2-20µl	
		✓	✓																				5-50µl	BrandTech Transferpette (single channel)
		✓	✓																				10-100µl	
		✓	✓											✓	✓	✓	✓						20-200µl	
		✓	✓											✓	✓	✓	✓						25-250µl	
				✓	✓	✓	✓	✓											✓	✓	✓		100-1000µl	
✓	✓									✓	✓	✓	✓										0.5-10µl	
				✓	✓																		2-20µl	
				✓	✓																		2.5-25µl	
				✓	✓																		5-50µl	BrandTech Transferpette (multichannel)
				✓	✓																		10-100µl	
				✓	✓																		20-200µl	
				✓	✓																		30-300µl	
✓	✓									✓	✓	✓	✓										0.5-10µl	
		✓	✓																				5-100µl	
		✓	✓											✓	✓	✓	✓						2-200µl	Eppendorf Electronic
				✓	✓	✓	✓	✓											✓	✓	✓		50-1000µl	



# Pipet Tip Compatibility

0.5-10µl	0.5-10µl	20-20010	1-200µl	21-0010	20-0200	1-200µl	21-0200	100-1000µl	20-1000	100-1000µl	21-1000	100-1000µl	20-1250	100-1000µl	20-1251	0.5-10µl	21-1250	0.5-10µl	22-0010	0.5-10µl	22-0011	0.5-10µl	23-0010	1-200µl	23-0011	1-200µl	22-0200	1-200µl	22-0201	1-200µl	23-0200	100-1000µl	23-0201	100-1000µl	22-1000	100-1000µl	22-1001	100-1000µl	23-1000	100-1000µl	23-1001	Volume/Biologix Tip	Volume Range	Pipettor Brand/Make														
✓	✓																																										0.1-10µl															
		✓	✓																																										2-20µl													
		✓	✓																																											10-100µl												
		✓	✓																					✓	✓	✓	✓																		20-200µl	Eppendorf Reference (single channel) & (multichannel)												
				✓	✓	✓	✓	✓	✓														✓	✓	✓	✓																			100-1000µl													
✓	✓																✓	✓	✓	✓																										0.1-2.5µl												
✓	✓																✓	✓	✓	✓																											0.5-10µl											
		✓	✓																																													2-20µl										
		✓	✓																																													10-100µl	Eppendorf Research (single channel)									
		✓	✓																																														25-200µl									
				✓	✓	✓	✓	✓																✓	✓	✓	✓																					100-1000µl										
✓	✓																✓	✓	✓	✓																														0.5-10µl								
																																																			5-100µl							
		✓	✓																																															50-100µl								
		✓	✓																																															20-300µl	Eppendorf Research (multichannel)							
		✓	✓																																															30-300µl								
		✓	✓	✓	✓																																														50-1200µl							
✓	✓																✓	✓	✓	✓																																0.2-10µl						
✓	✓																✓	✓	✓	✓																																0.5-10µl						
		✓	✓																																																	5-40µl	Thermo Scientific Finnpiptette (single channel)					
		✓	✓																																																40-200µl							
				✓	✓	✓	✓	✓																																													100-1000µl					
				✓	✓	✓	✓	✓																																														200-1000µl				
		✓	✓																																																		5-50µl	Thermo Scientific Finnpiptette (multichannel)				
		✓	✓																																																	50-300µl						
✓	✓																✓	✓	✓	✓																																		P2, 0.1-2µl				
✓	✓																✓	✓	✓	✓																																			P10, 0.2-10µl			
		✓	✓																																																					P20, 2-20µl	Gilson Pipetman (single channel)	
		✓	✓																																																				P100, 10-100µl			
		✓	✓																																																					P200, 20-200µl		
				✓	✓	✓	✓	✓																																																	P1000, 100-1000µl	
		✓	✓																																																						P200-M8, 20-200µl (m/ch)	

# Pipet Tip Compatibility

Pipettor Brand/Make		Volume Range	0.5-10µl	0.5-10µl	1-200µl	1-200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	0.5-10µl	0.5-10µl	0.5-10µl	0.5-10µl	1-200µl	1-200µl	1-200µl	1-200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl							
Volume/Biologix Tip			20-0010	21-0010	20-0200	21-0200	20-1000	21-1000	20-1250	21-1250	22-0010	22-0011	23-0010	23-0011	22-0200	22-0201	23-0200	23-0201	23-0201	23-0201	100-1000µl	100-1000µl	22-1000	22-1001	23-1000	23-1001	23-1001	23-1001	23-1001	23-1001	23-1001							
Rainin EDP3 (multichannel)	EP-M8-10, 0.5-10µl		✓	✓									✓	✓	✓	✓																						
	EP-M8-250, 10-250µl				✓	✓											✓	✓	✓	✓																		
Rainin Latch Mode	R2, 0.1-2µl		✓	✓									✓	✓	✓	✓																						
	R10, 0.2-10µl		✓	✓									✓	✓	✓	✓																						
	R20, 2-20µl				✓	✓																																
	R100, 10-100µl				✓	✓																																
	R200, 20-200µl				✓	✓												✓	✓	✓	✓																	
	R1000, 100-1000µl						✓	✓	✓	✓	✓																		✓	✓	✓							
Socorex Calibra/Accura (single channel)	0.1-2µl		✓	✓									✓	✓	✓	✓																						
	0.2-2µl		✓	✓									✓	✓	✓	✓																						
	0.5-10µl		✓	✓									✓	✓	✓	✓																						
	1-10µl		✓	✓									✓	✓	✓	✓																						
	2-20µl				✓	✓																																
	10-100µl				✓	✓																																
	20-200µl				✓	✓												✓	✓	✓	✓																	
100-1000µl						✓	✓	✓	✓	✓																		✓	✓	✓								
Socorex Calibra/Accura (multichannel)	1-10µl		✓	✓									✓	✓	✓	✓																						
	5-50µl				✓	✓											✓	✓	✓	✓																		
	10-100µl				✓	✓											✓	✓	✓	✓																		
	20-200µl				✓	✓											✓	✓	✓	✓																		
	50-200µl				✓	✓											✓	✓	✓	✓																		
VWR and VWR Signature (single channel)	0.1-2µl		✓	✓									✓	✓	✓	✓																						
	0.5-10µl		✓	✓									✓	✓	✓	✓																						
	2-20µl				✓	✓																																
	5-50µl				✓	✓																																
	10-100µl				✓	✓																																
	20-200µl				✓	✓											✓	✓	✓	✓																		
	50-250µl				✓	✓											✓	✓	✓	✓																		
100-1000µl						✓	✓	✓	✓	✓																		✓	✓	✓								
VWR and VWR Signature (multichannel)	1-10µl		✓	✓									✓	✓	✓	✓																						
	5-50µl				✓	✓																																
	20-200µl				✓	✓											✓	✓	✓	✓																		
	50-300µl				✓	✓											✓	✓	✓	✓																		

# Pipet Tip Compatibility

Pipette Brand		Eppendorf Research Plus					Gilson Pipetman(PG)					Gilson Pipetman(PN)					Gilson Pipetman(PL)									
Model / Range		0.1-2.5ul	0.5-10ul	2-20ul	10-100ul	20-200ul	100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	
Volume	Product No.	Universal Fit Filter Tips																								
0.5-10ul	23-0010S	✓	✓					✓	✓					✓	✓					✓	✓					
	23-0011S	✓	✓					✓	✓					✓	✓					✓	✓					
1-20ul	23-0020S			✓					✓							✓						✓				
	23-0021S			✓					✓							✓						✓				
1-50ul	23-0050S				✓					✓	✓					✓	✓					✓	✓			
	23-0051S				✓					✓	✓					✓	✓					✓	✓			
1-100ul	23-0100S				✓					✓	✓					✓	✓					✓	✓			
	23-0101S				✓					✓	✓					✓	✓					✓	✓			
1-300ul	23-0200S					✓					✓						✓						✓	✓		
	23-0201S					✓					✓						✓						✓	✓		
100-1000ul	23-1000S						✓					✓					✓						✓	✓		
	23-1001S						✓					✓					✓						✓	✓		

Pipette Brand		Thermo Finnpiquette (F1)								Thermo Finnpiquette (F2)								Thermo Finnpiquette (F3)															
Model / Range		0.2-2ul	1-10ul	1-10ul(µTip)	2-20ul	2-20ul(µTip)	5-50ul	5-50ul(µTip)	10-100ul	20-200ul	30-300ul	100-1000ul	0.2-2ul(µTip)	1-10ul	1-10ul(µTip)	2-20ul	2-20ul(µTip)	5-50ul	5-50ul(µTip)	10-100ul	20-200ul	100-1000ul	0.2-2ul	1-10ul	1-10ul(µTip)	2-20ul	2-20ul(µTip)	5-50ul	5-50ul(µTip)	10-100ul	20-200ul	100-1000ul	
Volume	Product No.	Universal Fit Filter Tips																															
0.5-10ul	23-0010S	✓	✓											✓										✓	✓								
	23-0011S	✓	✓											✓										✓	✓								
1-20ul	23-0020S				✓										✓											✓							
	23-0021S				✓										✓											✓							
1-50ul	23-0050S					✓			✓	✓						✓	✓									✓	✓						
	23-0051S					✓			✓	✓						✓	✓									✓	✓						
1-100ul	23-0100S					✓			✓	✓						✓	✓									✓	✓						
	23-0101S					✓			✓	✓						✓	✓									✓	✓						
1-300ul	23-0200S										✓						✓										✓	✓					
	23-0201S										✓						✓										✓	✓					
100-1000ul	23-1000S										✓							✓									✓	✓					
	23-1001S										✓							✓									✓	✓					

Pipette Brand		BIOHIT Proline								BIOHIT Proline (Plus)								BRAND																	
Model / Range		0.1-2.5ul	0.5-10ul	2-20ul	5-50ul	10-100ul	20-200ul	50-300ul	100-1000ul	0.1-3ul	0.5-10ul	2-20ul	5-50ul	10-100ul	20-200ul	30-300ul	100-1000ul	0.1-1ul	0.1-2.5ul	0.5-10ul	2-20ul	5-50ul	10-100ul	20-200ul	30-300ul	100-1000ul	0.1-1ul	0.1-2.5ul	0.5-10ul	2-20ul	5-50ul	10-100ul	20-200ul	30-300ul	100-1000ul
Volume	Product No.	Universal Fit Filter Tips																																	
0.5-10ul	23-0010S	✓	✓								✓	✓							✓	✓							✓	✓							
	23-0011S	✓	✓								✓	✓								✓	✓							✓	✓						
1-20ul	23-0020S			✓									✓															✓							
	23-0021S			✓									✓															✓							
1-50ul	23-0050S				✓	✓	✓	✓					✓	✓	✓	✓											✓	✓							
	23-0051S				✓	✓	✓	✓					✓	✓	✓	✓											✓	✓							
1-100ul	23-0100S					✓	✓	✓						✓	✓	✓											✓	✓							
	23-0101S					✓	✓	✓						✓	✓	✓											✓	✓							
1-300ul	23-0200S						✓	✓							✓	✓											✓	✓							
	23-0201S						✓	✓							✓	✓											✓	✓							
100-1000ul	23-1000S								✓								✓											✓	✓						
	23-1001S								✓								✓											✓	✓						