

Cryogenic

-150°C Chest Freezer

231 L



MDF-C2156VAN-PE

The most uniform storage temperatures for cryopreservation solutions

The MDF-C2156VAN Cryogenic Freezer offers unparalleled safety, reliability and uniform temperatures well below -135°C for ideal long-term cryopreservation solutions.

Optimum Sample Storage


The use of patented space saving VIP PLUS panels within the freezer cabinet provides a reduced wall thickness for maximum interior volume. This results in a capacity of up to 150 2" boxes.

Reliable Technologies

Compressors that are specifically designed for ultra-low temperature applications are employed in the proven cascade refrigeration system ensuring the highest levels of performance and reliability.


Ease of Use & Intelligent Security

All alarm functions, self-diagnostic notifications and a graphical display of temperature performance over time are available in the specially designed LCD control panel.



Uniform Sample Storage

The uniform and stable temperatures provided are ideal for storage of cells, tissues and other important biological samples.



Dry Storage

The highest levels of protection for valuable samples, without the risks of cross contamination and vertical temperature gradients.



Easy Monitoring

The blue display provides a clear view of the temperature and gives a notification in the case of abnormalities.

