



Well Organized

The new Eppendorf Tube and Cuvette Racks

Optimize your workflows through user-friendly handling of clearly arranged samples. Almost all laboratory protocols require efficient and reliable processing, transport, and short term storage of sample vessels in benchtop racks. The new Eppendorf Tube and Cuvette Racks combine optimized functionality and high robustness with an attractive design. Six different formats offer the optimal solution and perfect fit for all tubes and cuvettes commonly used in laboratories. Compact dimensions require minimal space.

- > All racks are stackable to save lab space when not in use
- > Made from high quality polypropylene (PP):
 - > Ensuring dimensional stability across a broad temperature range (-86 ° C to 121 ° C)
 - > Providing chemical resistance
- > The racks are laboratory washer-safe for fast and easy cleaning
- > The racks are autoclavable (121 ° C/20 min) and thus suitable for workflows requiring a sterile environment, e.g. cell handling applications



> Non-slip silicone feet and a tilt-proof footprint design ensure safe handling.




> A wavelike design enables you to clearly identify your samples, even when tubes are placed vertically in a column.



> Racks are compatible with all commonly used tubes and cuvettes.
> The 1.5 - 50 mL Tube Racks accommodate different tube formats on the same rack.



> All racks are stackable to save lab space when not in use.



> Consecutive numbering of rack positions supports clear allocation and error-free handling.



> The anti-rotation structure of the Cryogenic Tube Rack fits all commonly used cryotubes and enables one-hand operation.

Ordering information

Description

Eppendorf Tube Rack for 0.5 mL tubes, 48 wells (3 rows of 16 wells each), 2 pcs.
Eppendorf Tube Rack for 1.5/2.0 mL tubes, 36 wells (3 rows of 12 wells each), 2 pcs.
Eppendorf Tube Rack for 5.0/15 mL tubes, 12 wells (2 rows of 6 wells each), 2 pcs.
Eppendorf Tube Rack for 5.0/15/50 mL tubes, 12 wells (6 x 5.0/15 mL and 6 x 50 mL), 2 pcs.
Eppendorf Tube Rack for cryogenic tubes, 36 wells (3 rows of 12 wells each), 2 pcs.
Eppendorf Cuvette Rack, 30 wells (3 rows of 10 wells each), 2 pcs.

Order no.

0030 119.800
 0030 119.819
 0030 119.827
 0030 119.835
 0030 119.843
 0030 119.851

000 «Диаэм»

Москва
 ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
 +7 (812) 372-6040
 spb@dia-m.ru

Новосибирск
 +7 (383) 328-0048
 nsk@dia-m.ru

Воронеж
 +7 (473) 232-4412
 vrn@dia-m.ru

Йошкар-Ола
 +7 (927) 880-3676
 nba@dia-m.ru

Красноярск
 +7 (923) 303-0152
 krsk@dia-m.ru

Казань
 +7 (843) 210-2080
 kazan@dia-m.ru

Ростов-на-Дону
 +7 (863) 303-5500
 rnd@dia-m.ru

Екатеринбург
 +7 (912) 658-7606
 ekb@dia-m.ru

Кемерово
 +7 (923) 158-6753
 kemerovo@dia-m.ru

Армения
 +7 (094) 01-0173
 armenia@dia-m.ru

